


L&L SKIN

MIO 2

Important Reminder:

1. Please charge to the saturation before use.
2. Please recharge the device if it's not used for a long time (about 3 months).



Please charge the device before use:

It takes hours to full charge the device.

You can use the massager for 10 houses in blue mode or 30 min in red mode when fully charged.

The red-light blinking indicates low battery-please place the device into base for charging.

How to use the MIO2 Massager:

For the best result move the massager following your face muscles location (as shown on the picture).



The surface #1 of the device (as shown on the picture) is intended for Gua Sha massage. Gusa Sha is very popular in China and Japan-it allows you to improve blood circulation and realise discomfort. You can use for any part of your body.

The surface #2 is good for muscles around your cheeks, forehead and neck area.

The surface #3 is intended for muscles around your eyes, chin and lips.

1. Ensure that your skin is clean, and you don't have any make up on.
 2. Apply proper cream, lotion, serum for your type of skin
 3. Press on/off button for 1 second to switch on the massager.
 4. Briefly press the on/off button to switch between heated and soothing mode to achieve needed effect.
- Soothing mode (Blue light flash indication): cools and soothes your skin with a low vibration – 142 Hertz.
 - Heated mode (Red light flash indication): helps to rejuvenate your skin with higher vibration-200 Hertz and 45°C / 113°F heat.



5. The shapes of the device provide micro-currents and are smoothly-curved to enable easy application at any part of your face and body. Feel free to use any of the three devices surfaces, choosing the one that better fits you.
6. After the use, switch off the device by pressing the button for 1 second, wipe it out with a paper towel or wash it and dry in a well-ventilated place.

FAQ

For how long should I use the device every day to see the result?

We recommend using the device for 5-10 minutes every day after morning and evening washing. The best result you will see in 60 days of everyday use.

Can I apply a cream or serum on the skin before the procedure?

Yes, it is recommended to use cosmetic solution which fits your type of skin. It is also preferable to use the massager after moisturizing, collagen or vitamin masks-the device helps nutrients to get deeper into the skin.

Can I use the device every day?

Yes, the vibration intensity and temperature adjustment of this device are intelligently designed by a professional team. It is suggested to use multi times a day.

The product suddenly stopped working or doesn't switch on. What should I do?

- Check if the device was properly charged. Always recharge it before use if it was not in use for a long time.
- Check if the charging base connected to USB cable and the cable plugged into electricity source.
- The red-light blinking indicates low battery-please place the device into base for charging.
- Try to switch off and on the device once again, pressing the off/on button for one second.

Disclaimer

Not a medical device, should not be used to treat any specific diseases. Avoid using the device in case of any kind of skin pathology or disease, skin damage, wounds or scratches. In case if you feel any discomfort while using the massager you should stop using it immediately. Don't use the massager while charging.